

NEU


Implantate für die Hand

APTUS Hand-System 1.2 – 2.3

1.5 TriLock Implantate

Längere gerade und T-Platte

Längere TriLock Schrauben



2.0 TriLock Implantate

Rotationsplatte

Längere T-Platten




NEU

Bestellinformation

1.5 TriLock Platten

Material: Titanium (ASTM F67)
Plattendicke: 0.8 mm




A-4350.08

A-4350.17

Art. Nr.	Beschreibung	Löcher	Stk./Pkg.
A-4350.08	gerade	8	1
A-4350.17	T	10 (3/7)	1

2.0 TriLock Platten

Material: Titanium (ASTM F67)
Plattendicke: 1.3 mm



A-4655.11

A-4655.17

A-4655.24

Art. Nr.	Beschreibung	Löcher	Stk./Pkg.
A-4655.11	T	8 (3/5)	1
A-4655.17	T	10 (3/7)	1
A-4655.24	Rotation	6 (3/3)	1

1.5 TriLock Schrauben, HexaDrive 4

Material: Titanium (ASTM F136)



Länge	Art. Nr.	Stk./Pkg.	Art. Nr.	Stk./Pkg.
14 mm	A-5250.14/1	1	A-5250.14	5
16 mm	A-5250.16/1	1	A-5250.16	5
18 mm	A-5250.18/1	1	A-5250.18	5
20 mm	A-5250.20/1	1	A-5250.20	5